


Voltage
Signal 1/p
generator
with step
up transformer




Circuit Diagram






Ferro-electric crystal.

Sweyer Tower:



Expected graphs:



PE Hysterisis of Ferroelectric Crystal

Aim - To study the hysteresis loop of the given ferroelectric material and to observe the variation of parameters like coercive field, remnant polarization and saturated polarization.

Apparatus Required - Ferroelectric Sample, Sample Holder, Polarization vs electric field set-up, DSO.

Formula Used -

$$\textcircled{a} \quad \text{Coercive Field } (E_c) = \frac{K_I \times S_x}{t}, \quad K_I \text{ is intercept of PE loop}$$

S_x - Volts/division of graph

t - thickness of sample

$$\textcircled{b} \quad \text{Remnant Polarization } (P_r) = \frac{C \times Y_{\text{intercept of PE loop}} \times Y_{\text{volts/div}}}{{\text{Area of Sample}}}$$


C - Capacitance of standard capacitor.

$$\textcircled{c} \quad \text{Saturated Polarization } (P_s) = \frac{C \times Y_{\text{max of PE loop}} \times Y_{\text{volts/div}}}{{\text{Area of Sample}}}$$

Precautions - \textcircled{a} Handle apparatus with care, as we are dealing with High voltage.

\textcircled{b} Do not disturb the sample throughout the experiment.


1. $V = 503 V$




2. $V = 566 V$



3. $V = 617 V$



4. $V = 625 V$



Experiment:

33/3

Page No. 24

- (1) Make sure the sample is dipped in low viscosity oil (say Silicon oil) to prevent discharge due to High voltage.
- (2) Cover sample Holders with a lid while performing the experiment. This helps to have a stable temperature for the sample.

Observation: Sample dimensions : 5mm x 5mm

Thickness of sample : 0.5mm

Capacitance of standard capacitor, $C = 0.65 \text{ nF}$

Symbols used in observation table:

K_{IL}, K_{IR} - Left and Right X intercepts of PE loop respectively

K_I - Mean of K_{IL} & K_{IR}

y_{IU}, y_{ID} - Up and Down y-intercepts of PE loop respectively

y_I - Mean of y_{IU} and y_{ID}

y_{MU}, y_{MD} - Up & down y-coordinates of the maximum y-value in PE loop.

V - Applied voltage across the sample.

Experiment:

Date _____

Page No. 25

S.No.	V	$S_1(V_{de})$	$S_2(V_{de})$	X_{IL}	K_{IR}	$Y_{1,V}$	$Y_{1,D}$	Y_{MV}	Y_{MD}	K_I	Y_I	Y_M
1.	503	0.5	0.077	3.0	2.8	1.8	2.0	2.8	3.2	2.9	1.9	3.0
2.	566	0.5	0.074	3.2	3.2	2.2	2.4	3.4	3.8	3.2	2.3	3.6
3.	617	0.5	0.090	3.4	3.4	2.2	2.4	3.4	3.6	3.4	2.3	3.5
4.	625	0.5	0.124	3.4	3.6	1.6	1.8	2.4	2.6	3.5	1.7	2.5
5.	660	0.5	0.124	3.6	3.6	1.6	2.0	2.6	2.6	3.6	1.8	2.0
6.	700	1.0	0.124	1.8	2.0	1.8	2.0	2.8	3.0	1.9	1.9	2.9
7.	720	1.0	0.128	2.0	2.0	1.8	2.0	2.8	3.2	2.0	1.9	3.0
8.	760	1.0	0.128	2.0	2.0	2.0	2.2	3.0	3.4	2.0	2.1	3.2
9.	820	1.0	0.136	2.2	2.2	2.0	2.2	3.2	3.4	2.2	2.1	3.3
10.	840	1.0	0.140	2.2	2.2	2.0	2.2	3.2	3.4	2.2	2.1	3.3
11.	900	1.0	0.140	2.2	2.4	2.2	2.4	3.4	3.6	2.3	2.3	3.5
12.	940	1.0	0.152	2.4	2.4	2.0	2.2	3.4	3.6	2.4	2.1	3.5
13.	1000	2.0	0.220	1.2	1.4	1.6	1.8	2.6	2.8	1.3	1.7	2.7

Calculations: Area of Sample = $(5 \times 5) \text{ mm}^2 = 25 \times 10^{-6} \text{ m}^2$

S.No.	V (V)	$E_c(V/m)$	$P_t(k/m^2)$	$P_s(c/m^2)$
1.	503	2900	0.00363	0.00577
2.	566	3200	0.00443	0.00693
3.	617	3400	0.00538	0.00819
4.	625	3500	0.00548	0.00806
5.	660	3600	0.00580	0.00695
6.	700	3800	0.00613	0.00935
7.	720	4000	0.00632	0.00998
8.	760	4000	0.00699	0.01065
9.	820	4400	0.00743	0.01167
CLASSTIME				

1.	840	4400	0.00764	0.01201
2.	900	4600	0.00837	0.01274
3.	940	4800	0.00830	0.01383
4.	1000	5200	0.00972	0.01594

Result: Hysteresis loop of ferroelectric material is studied and variation of coercive field, remanent polarization and saturated polarization are plotted wrt applied sample voltage.

Discussion-X-Signal of PC setup represents voltage across sample while Y signal is proportional to charge on capacitor and Polarisation.

In general, polarization is proportional to field applied. But in case of ferro-electric there is remnant polarization, i.e. polarization remains in the direction of applied field even after this field is removed hence hysteresis is displayed.

When applied Voltage across sample is increasing, two things can be changed in ferroelectric material.

Swinger-Tower circuit:

It is used to analyse hysteresis loop for material which has low loss and high polarization. Sample is attached to series capacitance (C_s) where capacitance should be more than sample capacitance (C_{st}). C_s is reference capacitor that is known & independent.